


SOP-G3AC-01 · VERSION 1.0

StatPacks G3+ Airway Cell

Standard Operating Procedure for inspection, loadout, and field deployment. Adherence ensures operator proficiency, equipment readiness, and patient safety.



SECTION 1

Mission Brief & Intended Use

1.1 Primary Mission

Provide a standardized, protected, and rapidly deployable module for organizing BLS and ALS airway management equipment — reducing cognitive load and minimizing search time during high-stress respiratory emergencies.

1.2 Intended Use

Designed for use by **trained and certified medical personnel** — Paramedics, Flight Medics, EMTs, and Combat Medics — operating within their authorized scope of practice. **Not intended for use by untrained individuals.**

1.3 Operational Context

Integrates seamlessly into any StatPacks G3+ main pack (e.g., Backup, Breather). Serves as a dedicated "grab-and-go" airway kit, separating delicate airway tools from general trauma supplies to protect them from damage and contamination.

ALS Role

Organizes and protects equipment for advanced airway interventions:

- Direct laryngoscopy
- Endotracheal intubation
- Rapid Sequence Intubation (RSI)

BLS Role

Organizes and protects equipment for basic airway management:

- Oropharyngeal airways (OPAs)
- Nasopharyngeal airways (NPAs)
- Suction adjuncts

SECTION 2

Equipment & Component Familiarization

Components are grouped by operational purpose.

Equipment Protection & Structure

Semi-Rigid Foam-Molded Shell

Provides impact and crush resistance — critical for protecting laryngoscope blades, bulbs, fiber-optics, and pre-packaged sterile items if the main pack is dropped, compressed, or struck during operations.

Rapid Identification & Access

- **Green Exterior Color:** Universal EMS color code for airway equipment — enables immediate visual ID by any team member under duress.
- **High-Visibility Zipper Pulls:** Positive grip and rapid opening, even with gloved hands or in low-light environments.
- **Clear ID Window (Handle):** Insert custom labels (e.g., "ADULT ALS," "PEDIATRIC AIRWAY," "TRUCK 86") to eliminate confusion when multiple cells are stored together.

Interior Organization & Layout

- **Quad-Panel Clamshell Design:** Lays flat when opened, presenting equipment across four distinct panels. Supports muscle memory and allows simultaneous visualization of all tools.
- **Elastic Loops and Pockets:** Securely retains laryngoscope handles, blades (Macintosh/Miller), ET tubes, stylets, syringes, and airway adjuncts — preventing disorganization or damage during transport.

System Integration

Hook-and-Loop Rear Panel

Affixes the cell securely to the interior of a compatible StatPacks G3+ main pack. Prevents the module from shifting or falling out when the main pack is opened — ensuring it is always in its expected location.

Standard Operating Procedure: Deployment

- Scenario:** Unit dispatched to a 58-year-old female in respiratory arrest requiring immediate airway intervention. Partner is providing BVM ventilations while you prepare for endotracheal intubation.

1

Locate & Retrieve

- Open the primary medical pack (e.g., StatPacks G3+ Backup).
- Visually identify the Airway Cell by its **green color**.
- Grasp the handle and detach from the hook-and-loop mounting point. **Do not pull on the body of the cell.**

2

Stage for Intervention

- Place the cell on a stable, clean surface accessible to the airway provider — typically near the patient's head.
- Position so it can be opened fully without obstruction.

3

Deploy & Prepare Equipment

- Use zipper pulls to fully open the cell. Lay flat to expose all four internal panels.
- Select primary equipment: laryngoscope handle + blade (e.g., Mac 3), appropriately sized ET tube, stylet, and 10cc syringe.
- Assemble and functionally test all components. Check laryngoscope light and inflate/deflate ET tube cuff to check for leaks.
- Select backup equipment (smaller/larger ET tube, alternate blade, OPA).

4

Post-Intervention Recovery

- Properly discard all disposable materials.
- Segregate all reusable equipment (laryngoscope handle, blades) for decontamination per service protocol.
- Secure any remaining unused supplies within the cell.

5

Reconstitution

The cell **must be fully restocked and decontaminated** before returning to operational status. **Do not place a used or partially stocked cell back into the primary pack.**

SECTION 4

Critical Warnings, Limits & Safety

Scope of Practice Warning

This module holds tools for advanced medical procedures carrying a **high risk of patient harm — including hypoxia, trauma, and death — if performed improperly**. Such interventions must only be attempted by currently certified, credentialed, and authorized medical professionals acting within their local protocols and scope of practice.

Equipment, Not a Medical Device

The G3+ Airway Cell is an **organizational tool only**. It provides no therapeutic benefit on its own. Efficacy of any intervention depends entirely on the skill of the provider and the condition of the medical supplies stored within.

Physical Limitations

The semi-rigid shell provides significant protection against drops and compression but **is not indestructible**. It will not protect contents from extreme crushing forces (e.g., being run over by a vehicle) or submersion in water.

Cross-Contamination Hazard

Airway management generates aerosols and involves contact with bodily fluids. The cell's interior and exterior must be considered **contaminated after every use**. Follow service-specific protocols for decontamination using an approved virucidal/bactericidal agent to prevent pathogen transmission.

SECTION 5

Readiness, Inspection & Sustainment

5.1 Pre-Shift / Daily Inspection

Exterior:

- Inspect cell for cuts, tears, or shell damage.
- Verify zipper pulls are present and zippers function smoothly along their entire length.

Interior:

- Open cell and conduct full inventory against service's prescribed load list.
- Check for damaged packaging on sterile items.

Contents:

- Confirm laryngoscope functionality (light source).
- Verify all medications (if carried) and supplies are within expiration dates.
- Ensure all items are properly secured in designated loops/pockets.

5.2 Post-Mission Sustainment

Restock: Immediately after any operational use, take the cell out of service until a complete inventory and restock is performed. **Every single item used must be replaced.**

Decontaminate: Clean and decontaminate interior and exterior surfaces per bloodborne pathogen and infectious disease control protocols. Pay special attention to high-contact areas: handle and zippers.

5.3 Periodic Maintenance

On a **monthly or quarterly basis:**

- Fully empty the cell and clean the interior thoroughly.
- Inspect elastic loops for loss of elasticity.
- Inspect interior lining for signs of wear or degradation.