G3+ Clinician Pro

Field Manual & Standard Operating Procedure

STATPACKS G3+ Clinician Pro ALS backpack
Red G35001PRE or Blue G35001PBU, 30 liters



Section 1 Mission Brief & Intended Use	**Section 2** Component Familiarization	**Section 3** Deployment SOP	**Section 4** Warnings & Limits	**Section 5** Readiness & Sustainment

Prepared under the authority of Valerie Ellis Lavin, Founder and CEO, Retired U.S. Army First Sergeant

This field manual supplements, and does not replace, the manufacturer's instructions and labels. Follow them, and your agency or organization's procedures where they apply. Product names are trademarks of their respective owners.

SECTION 1

Mission Brief & Intended Use

1.1 Mission Objective

Carry a full ALS setup in a compact pack built around three Generation 3 Cells, with a shelf layout you set once and repeat.

1.2 Scope of Application

Covers inspecting, configuring, loading, carrying, working from and restocking the pack. Cells, kits, supplies and their use follow your agency's protocol and the training each medic holds.

The STATPACKS G3+ Clinician Pro is a compact 30-liter ALS backpack for medics, sized to carry three Generation 3 Cells on an adjustable track-and-shelf system, with a vertical design that keeps the straps off the ground.

1.3 OPERATIONAL CONTEXTS

EMS and fire crews

ALS medics who want three Cells and nothing extra, in a 30-liter backpack.

Water rescue and sports

Lifeguard, water rescue and sports medicine staff setting the shelves to fit three Cells.

Tactical medic teams

Medics working from an upright pack that keeps its straps elevated and off the ground.

SECTION 2

Component Familiarization

Know every component by name and by feel before you need it.

2.1 INTERIOR

Adjustable track-and-shelf system

Included. Move or reposition shelves to fit three Generation 3 Cells, which are not included.

Main compartment

21" x 12" x 8", sized around three Generation 3 Cells in a vertical layout.

2.2 OUTSIDE

Side pockets and top pocket

Side pockets 14" x 7" x 3" and a top pocket 12" x 6" x 3" for frequently used items.

2.3 BUILD AND CARRY

Vertical design and straps

Keeps the straps elevated and off the ground while you access contents on scene.

PVC-free tarpaulin

Non-toxic and urethane coated. 22" H x 18" W x 8" D, 30 liters, 4.8 lbs empty, Red or Blue.

SECTION 3

Standard Operating Procedure: Deployment

Scenario: During a station readiness check, a medic moves a shelf on the track-and-shelf system to fit the service's three Generation 3 Cells. She stages frequently used items in the side pockets, closes the pack and returns it to its compartment on the unit.

01

3.1 Prepare

1. **Inspect on arrival.** Confirm 1 G3+ Clinician Pro backpack and its adjustable track-and-shelf system, in the color you ordered, Red or Blue. Check the straps, side pockets, top pocket and tarpaulin. The pack ships empty, with no Cells or kits.
2. **Plan the loadout.** Build the contents from your agency or department loadout protocol around three Generation 3 Cells, which are not included. Give each Cell a shelf position and list the frequently used items for the side pockets and top pocket.
3. **Set the shelves.** Move or reposition shelves on the track-and-shelf system to fit the Cells and gear you carry. Set every pack on the service to the same layout, so each Cell sits where any medic's hands expect it.

02

3.2 Load and carry

4. **Load the pack.** Seat the Cells on their shelves in the 21" x 12" x 8" main compartment, heavier items lower as the shelf layout allows. Stage frequently used items in the side pockets and the top pocket.
5. **Close and check.** Close the main compartment and every pocket fully without forcing them over the contents. If a closure strains, reset a shelf or move an item, then close it again.
6. **Carry to scene.** Carry the Clinician Pro on its straps as a backpack. It weighs 4.8 lbs empty, so the loadout you add sets the carry weight; pack only what you can carry safely to the patient.

03

3.3 Work and recover

7. **Work from the pack.** Set the pack down at your side and open it. The vertical design keeps the straps elevated and off the ground while you access contents. Use what it carries only as your training and agency protocol direct.
8. **Unpack and restock.** After the call, remove used items and dispose of them as your protocol directs. Clean the pack as your agency directs, restock the Cells and pockets to the loadout list, check any dates printed on the packaging and restage it.

INSPECT	PLAN	SHELVE	LOAD	CARRY	RESTOCK

SECTION 4 · READ CAREFULLY

Critical Warnings, Limits & Safety

Use note: The Clinician Pro ships empty with its track-and-shelf system. Cells, kits, supplies and their use follow your agency's protocol and the training each medic holds.

4.1 HAZARDS

Load weight

Pack only the load you can carry safely for the distance and terrain, and keep heavy Cells seated on their shelves.

Overpacking

Do not force a pocket or compartment closed over its contents. Remove an item or reset the shelves instead.

Unseated shelves

A shelf that is not set on its track can shift under load. Confirm each shelf sits firmly before you load Cells.

4.2 LIMITS OF USE

Pack only

The pack carries what you load. Cells and the Intravenous, Medicine and Airway kits and supplies are not included.

When to escalate

If a strap, shelf or pocket fails during a shift, move the load to a serviceable pack and report it to your supervisor.

4.3 EQUIPMENT INTEGRITY

Shelves and track

Before each shift, check that each shelf moves, sets and holds on the track-and-shelf system.

Straps and tarpaulin

Check the straps and the urethane-coated tarpaulin for tears or wear, and retire the pack when they fail.

SECTION 5

Readiness, Inspection & Sustainment

5.1 Staging for Deployment

- Stage the loaded pack in the unit compartment or station spot your agency assigns, reachable without unloading other gear.
- Keep each Cell in the same shelf position in every pack on the service.
- Post a loadout list for the side pockets and top pocket and check it at each shift change.
- Issue packs by color, Red or Blue, if your agency assigns them by role or unit.

5.2 Periodic Inspection

- [ ] **Shelves**

  Confirm every shelf is set on its track and holds its Cell firmly.

- [ ] **Contents**

  Match each Cell and pocket to the loadout list and check any dates printed on the packaging.

- [ ] **Carry and shell**

  Check the straps, pockets and tarpaulin for fraying, cuts or residue.

5.3 SERVICE LIFE & REPLACEMENT

> Restock after every call and inspect the pack at each shift change. Keep it in service while the straps, shelves and tarpaulin are sound, and retire it when any of them fail.

This field manual supplements, and does not replace, the manufacturer's instructions and labels. Follow them, and your agency or organization's procedures where they apply. Product names are trademarks of their respective owners.